




**FIG.1**

**FIG.2**

**FIG.3**



**FIG. 6**

TABLE

PARAMETER	POSSIBLE RANGE	TYPE
SPHERE	-25 D TO +10 D	INPUT
CYLINDER	-10 D TO +10 D	INPUT
AXIS	0° TO 180°	INPUT
K-READING	33 D TO 53 D	INPUT
OPTICAL ZONE	4.5mm TO 8.5mm	INPUT
FLAP THICKNESS	0 μ m TO 300 μ m	INPUT
PACHYMETRY	300 μ m TO 700 μ m	INPUT
PUPILLOMETRY	4mm TO 11mm	INPUT
MAX. ABLATION	UP TO 250 μ m	OUTPUT
NUMBER OF PULSES	UP TO 9000	OUTPUT

FIG. 19

Smith, John OD: Date: _____		Treatment	
Less (μm): 180 Optical Zone (mm): 6.5 Pachymetry (μm): 523			
Information Rotational Eye Tracker Information Recommendation Information Restore Defaults			
<input type="radio"/> B Zyoptix B	Zyoptix B NORUS (6mm): 0.32 μm Pachymetry: 523 μm Acoustic Factor: 0.92+(Pac + 0) Subjective Refraction: -2.0/0.5/10° Zywave PPR: 2.47/0.63/11° Zywave Pupil Diameter: 6mm Treatment Sphere: -2.47 dpt Cylinder: 0.63 dpt Axis: 11° Number Of Pulses: 1744 Treatment Time: ??? Calculated Central Ablation: 45 μm Max Ablation: 47 μm Remaining Stroma: 313 μm Optical Zone: 6.50 mm		
<input type="radio"/> C Zyoptix C	Zyoptix C Max Ablation: 50 μm Remaining Stroma: 303 μm Optical Zone: 6.50 mm Treatment Sphere: -2.47 dpt Cylinder: 0.63 dpt Axis: 11° Number Of Pulses: 2144 Treatment Time: ??? Calculated Central Ablation: 45 μm PreOp K-Reading: 44.1 dpt PreOp Conic Const: -0.17 Pachymetry: 523 μm Acoustic Factor: 0.92-(Pac + 0) Subjective Refraction: -2.0/0.5/10° Treatment Sphere: -2.0 dpt Cylinder: 0.5 dpt Axis: 10° PreOp K-Read: 44.1 dpt PreOp Conic: -0.17 Desired: -0.56 Number Of Pulses: 2052 Treatment Time: ???		
<input type="radio"/> D Zyoptix D	Zyoptix D Max Ablation: 46 μm Remaining Stroma: 314 μm Optical Zone: 6.50 mm Treatment Sphere: -2.0 dpt Cylinder: 0.5 dpt Axis: 10° PreOp K-Read: 44.3 dpt PreOp Conic: -0.17 Desired: -0.56 Number Of Pulses: 1922 Treatment Time: ???		
<input type="button" value="Calculate"/>			
<input type="button" value="Zyoptix Export"/>			
<input type="button" value="Print"/>			

FIG. 7




FIG. 8

BAUSCH & LOMB OS 03/29/2001 Doe, John

Birthdate: 12/28/1988 Data Check Treatment Planner Version 0.02

Save Cancel

Pachymetry
 User Selected Pachymetry:
 Orbscan Pachymetry: 523 μ m
 Orbscan Exam Date: 12/28/2002 12:33:45 PM
 Remaining Pachymetry Under Flap: ???
 Flap Thickness: 160 μ m

Zywave Exam Information
 Zywave Exam Date: 12/28/2002 12:33:45 PM
 Zywave PPR (6mm): 2.47 0.63/11
 Zywave RMS (6mm): ???
 Rotational Eye Tracker Used

Zyoptix A Treatment 03/21/2003 2:32:22 PM
 Sphere: -2.47 dpt
 Cylinder: -0.63 dpt
 Axis: 11
 Manifest Refraction: ???
 Number Of Pulses: 2144
 Laser Frequency: 50 Hz
 Treatment Time: ???
 Optical Zone: 6.5 mm
 Treatment Zone: ???
 Max Ablation: 50 μ m
 Central Ablation: 45 μ m

315a

< | ? Demo Version

FIG.9A

BAUSCH & LOMB OD 01/23/2002 Doe, Jane

Birthdate: 12/28/1988 Data Check Treatment Planner Version 0.02

Save Cancel

Pachymetry
 User Selected Pachymetry:
 Orbscan Pachymetry: 523 μ m
 Orbscan Exam Date: 12/28/2002
 Remaining Pachymetry Under Flap: ???
 Flap Thickness: 160 μ m

Zyoptix C Treatment 03/21/2003 2:32:22 PM
 Sphere: -2.47 dpt
 Cylinder: -0.63 dpt
 Axis: 11
 PreOp K-Read: 44.1 dpt
 Q Pre: -0.17
 Q Post: ???
 Manifest Refraction: ???
 Number Of Pulses: 2144
 Laser Frequency: 50 Hz
 Treatment Time: ???
 Optical Zone: 6.5 mm
 Treatment Zone: ???
 Max Ablation: 50 μ m

315b

< | ? Demo Version

FIG.9B

Defaults

Restore Defaults

Set

Cancel

Limits

Zyoptix A

Min Remaining Stroma (μm):

Max Number of Pulses:

Max Treatment Time:

Max Ablation (μm):

Min Optical Zone (mm):

Zyoptix B

Min Remaining Stroma (μm):

Max Number of Pulses:

Max Treatment Time:

Max Ablation (μm):

Min Optical Zone (mm):

Zyoptix C

Min Remaining Stroma (μm):

Max Number of Pulses:

Max Treatment Time:

Max Ablation (μm):

Min Optical Zone (mm):

Max PostOp K-Reading (dpt):

Min PostOp K-Reading (dpt):

Zyoptix D

Min Remaining Stroma (μm):

Max Number of Pulses:

Max Treatment Time:

Max Ablation (μm):

Min Optical Zone (mm):

Max PostOp K-Reading (dpt):

Min PostOp K-Reading (dpt):

Preferences

◁ || ?

316

FIG.10




FIG.11A




FIG.11B

BAUSCH & LOMB

Preferences

Settings

General

Ote File Path: Language:

Ate File Path: Laser Frequency:

TLS File Path:

Clinic

Clinic Name:

Software Location:

Laser Number: Surgeon Name: Technician Name:

Zyoptix Treatment Planner

1202a

FIG.12A

BAUSCH & LOMB

Preferences

Defaults

Flap Thickness (μm):

Zyoptix Tissue Saving:

Optical Zone (mm)

☐ Use Default Value Default Value:
☒ Calculate from Pupillometry Add Value to Pupillometry:
Stay within this range. Min: Max:

Percent of Baseline Refraction

Procedure	Refraction Type	Sphere	Cylinder	
Zyoptix Personalized	Hyperopia	100	100	<input type="button" value="v"/>
Zyoptix Personalized	Myopia	100	100	<input type="button" value="v"/>
Zyoptix Personalized	Hyperopia	100	100	<input type="button" value="v"/>
Zyoptix Personalized	Myopia	100	100	<input type="button" value="v"/>

Zyoptix Treatment Planner

1202b

FIG.12B

BAUSCH & LOMB

1302

1318

1310

1304

1306

1308

1312

1316

1304

Patient Search

Last Name

First Name

OTE Files

Smith, John

OD

03/18/2003 11:43:22 AM

DOB (mm/dd/yyyy):

Gender:

Sub Refraction:

File Name:

Comment:

ATE Files

Smith, John

OD

09/17/2003 16:34:25

Patient ID:

DOB (mm/dd/yyyy):

Gender:

Sub. Ref.:

PPR:

Wavefront Dia.(mm):

File Name:

Comment:

Zyoptix Treatment Planner

FIG.13

BAUSCH & LOMB

1402

Continue

Restore

Cancel

Patient Information

Patient ID

000011

Last Name

Smith

First Name

John

DOB (mm/dd/yyyy)

05/27/1967

Age

36

Gender

☒ Male ☐ Female

Eye Information

Eye

☒ OD ☐ OS

☐ Retreatment

Subj. Sphere

-2.50

PPR Sphere

-2.39

Subj. Cylinder

-1.25

PPR Cylinder

-1.32

Subj. Axis

170

PPR Axis

163

Pachymetry (μm)

616

Pupillometry (mm)

6.4

Zyoptix Treatment Planner

FIG.14

BAUSCH & LOMB

Smith, John OD

Options

Personalized

Tissue Saving

Maps

Calculate

1302

Treatment Overview

Subj. Refraction: -2.50/-1.25/170

PPR (3.5mm): -2.39/-1.32/163

Wavefront Dia.: 7.01 mm

Pupillometry: 6.4 mm

Pachymetry: 616 μ m

INPUT VARIABLES

OPTICAL ZONE: 8.0

FLAP THICKNESS: 160

Pre-Op (6mm)

HO-RMS: 0.48 μ m

Z400: -0.38 μ m

Q (Conic): -0.21

K-Reading: 43.5 D

1504

1508

	Personalized	Tissue Savings
Baseline Ref.	-2.39/-1.32/163	-2.50/-1.25/170
Treatment Ref.	-2.39/-1.32/163	-2.50/-1.25/170
% of Baseline:	100%/100%	100%/100%
Optical Zone:	6 mm	6 mm
Flap Thickness:	160 μ m	160 μ m
Maximum Abl:	54 μ m	52 μ m
Central Abl:	54 μ m	52 μ m
Res. Stroma:	395 μ m Est.	397 μ m Est.
Treat. Zone:	8.3 x 8.2 mm	9.3 x 8.4 mm
Pulses (Time):	1858 (19 sec)	1908 (19 sec)
	Export TLB	Export TLB

Messages:

1508

Zyoptix Treatment Planner

FIG.15

BAUSCH & LOMB

Smith, John OD

Options

Calculate

Zyoptix Export

1615

1604, 1606, 1608

1602

Zyoptix Tissue Saving Treatment

Pre-Op K-Read (D): 43.5 Default

Pupillometry:

Pachymetry: 616 μ m

Treatment

Subj. Refraction: -2.50 / -1.25 / 170

Treatment Refraction: -2.50 / -1.25 / 170

Percent Subj. Ref.: 100 / 100

Optical Zone (mm): 8.0

Flap Thickness (μ m): 160

Max Ablation: 52 μ m

Residual Stroma: 397 μ m Est.

Treatment Zone: 9.3 x 8.4 mm

Pulses (Time): 1908 (19 sec)

Calc.

1612

1610

Ablation Profile

48

42.5

38.9

31.4

25.8

20.3

14.8

9.2

3.7

Messages:

1612

Zyoptix Treatment Planner

FIG.16




FIG. 17




FIG. 18A

BAUSCH & LOMB Smith, John OD

Data Check

1804

1802b

Save
Print
Save and Print
Treatment Overview
Select Patient
Finished

Laser Frequency: 100 Hz
Treatment Type: Myopia
PPR Type: Myopia
Subjective Refraction: -2.50/-1.25/170
PPR (3.5mm): -2.39/-1.32/163
Treatment Refraction: -2.39/-1.32/163
Percent of PPR: 100%/100%

Topo Thinnest Pachy: ---
Entered Pachymetry: ---
Flap Thickness: 160 μ m
Est. Residual Stroma: 395 μ m
Estimated from Post-Op Pachymetry Map and Flap Thickness.

Optical Zone: 6 mm
Max. Ablation: 54 μ m
Central Ablation: 54 μ m
Treatment Zone: 8.3 x 8.2 mm
Number Pulses (Time): 1858 (19sec)
Pulses for 2mm/1mm: 1317/539
Num. Treatment Phases: 2

Iris Recognition: Available
Ave Pupil shift X/Y: 28 μ m / 80 μ m
Ave Angle dil. vs undil.: -0.35°

TLB file:

Messages
Recommended Centration is the Pupil Center.

Entered Pupillometry:
Wave Pupillometry: 6.4 mm
Wavefront Diameter: 7.01 mm
HO-RMS (6mm): 0.49 μ m
Z400 (6mm): -0.36 μ m

Zyoptix Treatment Planner

FIG.18B

BAUSCH & LOMB Smith, John OD

Data Check

1802c

Save
Print
Save and Print
Treatment Overview
Select Patient
Finished

Pulses for 2mm/1mm: 1317/539
Num. Treatment Phases: 2

TLB file:

Messages
Recommended Centration is the Pupil Center

Print Comments

Visual Acuity SC:
Visual Acuity CC:
Planned Ring Size:

Entered Pupillometry:
Wave Pupillometry: 6.4 mm
Wavefront Diameter: 7.01 mm
HO-RMS (6mm): 0.49 μ m
Z400 (6mm): -0.36 μ m

Pre-Op Q (Conic,6mm): -0.21
Pre-op K-Reading(6mm): 45.5 D
White to White: ---
Simks: 46 D & 44.6 D
Topo Astig.: -1.4 D @ 168°
Angle Kappa: ---
Pup. Center vd Apex: ---

Zyoptix Treatment Planner

FIG.18C